

The Baxter logo is centered in the upper half of the slide. It consists of the word "Baxter" in a bold, italicized, blue sans-serif font. The background of the slide is a complex geometric pattern of overlapping diagonal bands in various shades of blue, gold, and light grey, creating a sense of depth and movement.


***Baxter***

IQE Reimagined


# Why Update IQE?

Creating design consistency across applications is extremely important to improving Baxter's brand identity as well as overall User Experience. As new doctors, nurses and clinicians come into the medical field they will be expecting clean, sleek and easy to use designs. IQ Enterprise is a very powerful application but unfortunately its visual design is rather dated. We took the new Baxter branding as well as the design concepts behind "Fugu" and did a "reimagining" of what IQE could look like with just a few changes. We took a few IQE screens and applied Fugu styling to see how they could be improved.

## Current IQE



## Updated IQE



# Update Process

---

The design change process was fairly simple. We looked at what made Fugu a visually appealing design and tried to transfer those elements to IQE.

- We added a more subtle gray background
- Used the same font styles
- Added newer icons where it made sense
- Experimented with different color combinations based on the new marketing guidelines
- Moved some page elements to help create a more consistent hierarchy
- Updated button and text link styles
- Recreated table and sort styling

Even though we made several visual updates and some interaction changes we did not change the core functionality of IQE.

- The navigation remains intact
- Data display is almost identical
- Filtering is untouched

# Home/Landing Page

The screenshot shows the home page of the IQ Enterprise Connectivity Suite. At the top left is the Baxter logo. To the right of the logo is a help icon (question mark) and the text "IQ Enterprise CONNECTIVITY SUITE". Further right is a user profile icon with the initials "AB". Below the header is a dark blue banner with the text "IQ Enterprise Connectivity Suite" and a sub-header: "IQ Enterprise Connectivity Suite is a centralized platform for Baxter applications and devices that help inform, improve and deliver hospital therapies." Below the banner is a search bar with the placeholder text "Search" and a magnifying glass icon. Under the search bar are two tabs: "ALL" (which is selected) and "INFUSION PUMP". Below the tabs is a section titled "INFUSION PUMP APPLICATIONS" with a light blue background. This section contains four application cards: 1. "Configuration Management" with a sub-header "Configure device distributions" and an icon of a central 'Rx' with four arrows pointing outwards. 2. "Device Management" with a sub-header "Device registration and management" and an icon of a network diagram with a central node. 3. "Reports" with a sub-header "CQI Reports and analytics" and an icon of a document with a checkmark and a clock. 4. "Enterprise Hierarchy" with a sub-header "Manage Enterprise Hierarchy" and an icon of a tree diagram. At the bottom left of the page is the copyright notice "© 2022 Baxter". At the bottom right are the version number "3.0.1" and two links: "Privacy Policy" and "Terms & Conditions".

Current IQE

# Home/Landing Page

**Baxter** IQ Enterprise CONNECTIVITY SUITE JD

## IQ Enterprise Connectivity Suite

IQ Enterprise Connectivity Suite is a centralized platform for Baxter applications and devices that help inform, improve and deliver hospital therapies.

### Applications

ALL INFUSION PUMPS Search

- Configuration Management**  
Configure device distributions
- Device Management**  
Device registration & Management
- Reports**  
CQI Reports and Analytics
- Enterprise Hierarchy**  
Manage Enterprise Hierarchy


© 2022 Baxter 3.0.1 | [Privacy Policy](#) | [Terms & Conditions](#)

## IQE Baxter Blue

Following the new brand guidelines we used the Baxter Blue (#154898) for all text links and buttons. While we used the navy blue (#001A72) for the header and titles.

- We added more rounded edges for the form fields.
- Reduced the size of the search bar
- Converted top tabs into a toggle
- Main Page Buttons
  - Removed the cyan
  - Updated some icons
  - Added subtle drop shadows

# Home/Landing Page




## IQE Navy Blue

We wanted to experiment and show the two branding blues in different roles.


For this design we used the navy blue (#001A72) for all text links and buttons. While we used the Baxter Blue (#154898) for the header and titles.

# Home/Landing Page: Navigation



Current IQE

# Home/Landing Page: Navigation



## IQE Baxter Blue

Updating the navigation was important. The “Infusion Pumps” dropdown was confusing when there were no other options.

Also, the buttons blended together as there were no separators.

We thought adding a black overlay would help differentiate the navigation from the rest of the.

Lastly, the text labels did not line up, for the update we lined them up to make them easier to read.



# Home/Landing Page: Navigation

The screenshot shows the top navigation bar of the IQ Enterprise Connectivity Suite. It features the Baxter logo on the left, a help icon, the product name 'IQ Enterprise CONNECTIVITY SUITE', and a user profile icon 'JD'. Below the navigation bar is a dark blue header with the title 'IQ Enterprise Connectivity Suite' and a brief description. The main content area has a light blue background with a search bar and a filter for 'ALL' and 'INFUSION PUMPS'. There are four main navigation cards: 'Configuration Management' (Configure device distributions), 'Device Management' (Device registration & Management), 'Reports' (CQI Reports and Analytics), and 'Enterprise Hierarchy' (Manage Enterprise Hierarchy). The footer contains copyright information and links to the Privacy Policy and Terms & Conditions.

## IQE Navy Blue

Updating the navigation was important. The “Infusion Pumps” dropdown was confusing when there were no other options.

Also, the buttons blended together as there were no separators.

We thought adding a black overlay would help differentiate the navigation from the rest of the.

Lastly, the text labels did not line up, for the update we lined them up to make them easier to read.

# Configuration Management

**Baxter** | IQ Enterprise CONNECTIVITY SUITE | AB

## Configuration Management

Configurations | Distributions

Configuration Type:

Search By:

<input type="checkbox"/>	Library Name	Device Model	Device Type	Library Version	Devices in Use	Imported By	Import Date
<input type="checkbox"/>	GatewayCQITestEvoLVP	Evo IQ	LVP	3	0	administrator	10-MAR-2022, 10:04
<input type="checkbox"/>	NovumLVPSYR_ROW2	Novum IQ	LVP	2	0	administrator	14-MAR-2022, 10:05
<input type="checkbox"/>	GatewayCQITestEVO	Evo IQ	LVP	2	0	administrator	10-MAR-2022, 10:03
<input type="checkbox"/>	GatewayCQITestEVOSYR	Evo IQ	SYR	1	0	administrator	10-MAR-2022, 10:04
<input type="checkbox"/>	GatewayCQI Test	Novum IQ	LVP	3	0	administrator	10-MAR-2022, 10:02
<input type="checkbox"/>	GatewayCQI Test	Novum IQ	LVP	5	0	administrator	10-MAR-2022, 10:03
<input type="checkbox"/>							10-JUN-2022

10 Items per page | 1 2 | 15 Total items |

Current IQE

# Configuration Management

**Baxter** IQ Enterprise CONNECTIVITY SUITE JD

## Configuration Management

Configurations Distributions

Select Configuration Type: Drug Library Search By: Search By Search

<input type="checkbox"/>	Library Name	Device Model	Device Type	Library Version	Devices in Use	Imported By	Import Date															
<input type="checkbox"/>	Very Long Library Name	Evo IQ	Syringe	1.2.3.4.5.6.7.8.9.10	100000	John Jacob Smith	MM/DD/YYYY															
<input type="checkbox"/>	Very Long Library Name	Evo IQ	Syringe	1.2.3.4.5.6.7.8.9.10	100000	John Jacob Smith	MM/DD/YYYY															
<table border="1"><thead><tr><th>Enterprise Unit Path</th><th># with Version</th><th>% with Version</th></tr></thead><tbody><tr><td>South West &gt; California &gt; Los Angeles &gt; LA Regional Hospital &gt; <b>Pediatrics</b></td><td>3</td><td>100%</td></tr><tr><td>South West &gt; California &gt; San Francisco &gt; SF Regional Hospital &gt; <b>Pediatrics</b></td><td>2</td><td>100%</td></tr><tr><td>South West &gt; California &gt; Sacramento &gt; Sacramento Regional Hospital &gt; <b>Pediatrics</b></td><td>2</td><td>100%</td></tr><tr><td>South West &gt; California &gt; San Diego &gt; San Diego Regional Hospital &gt; <b>Pediatrics</b></td><td>2</td><td>100%</td></tr></tbody></table>								Enterprise Unit Path	# with Version	% with Version	South West > California > Los Angeles > LA Regional Hospital > <b>Pediatrics</b>	3	100%	South West > California > San Francisco > SF Regional Hospital > <b>Pediatrics</b>	2	100%	South West > California > Sacramento > Sacramento Regional Hospital > <b>Pediatrics</b>	2	100%	South West > California > San Diego > San Diego Regional Hospital > <b>Pediatrics</b>	2	100%
Enterprise Unit Path	# with Version	% with Version																				
South West > California > Los Angeles > LA Regional Hospital > <b>Pediatrics</b>	3	100%																				
South West > California > San Francisco > SF Regional Hospital > <b>Pediatrics</b>	2	100%																				
South West > California > Sacramento > Sacramento Regional Hospital > <b>Pediatrics</b>	2	100%																				
South West > California > San Diego > San Diego Regional Hospital > <b>Pediatrics</b>	2	100%																				
<input type="checkbox"/>	Very Long Library Name	Evo IQ	Syringe	1.2.3.4.5.6.7.8.9.10	100000	John Jacob Smith	MM/DD/YYYY															

10 Items per page < 1 2 3 ... 8 9 10 > 400 Total items Next Step

© 2022 Baxter 3.0.1 | [Privacy Policy](#) | [Terms & Conditions](#)

## IQE Baxter Blue

The first thing to notice is the design contains an expandable table row, while what has been built does not. We thought this would give an opportunity to explore additional ideas.

The updated design follows the white background for the headline, gray body section, white pagination and gray footer.

Updates:

- We restructured the configuration and search areas
- There is vertical spacing between table rows, this makes it easier to scan and looks “cleaner”

# Configuration Management

**Baxter** IQ Enterprise CONNECTIVITY SUITE JD

## Configuration Management

Configurations Distributions

Select Configuration Type: Drug Library Search By: Search By Search

<input type="checkbox"/>	Library Name	Device Model	Device Type	Library Version	Devices in Use	Imported By	Import Date															
<input type="checkbox"/>	Very Long Library Name	Evo IQ	Syringe	1.2.3.4.5.6.7.8.9.10	100000	John Jacob Smith	MM/DD/YYYY															
<input type="checkbox"/>	Very Long Library Name	Evo IQ	Syringe	1.2.3.4.5.6.7.8.9.10	100000	John Jacob Smith	MM/DD/YYYY															
<table border="1"><thead><tr><th>Enterprise Unit Path</th><th># with Version</th><th>% with Version</th></tr></thead><tbody><tr><td>South West &gt; California &gt; Los Angeles &gt; LA Regional Hospital &gt; Pediatrics</td><td>3</td><td>100%</td></tr><tr><td>South West &gt; California &gt; San Francisco &gt; SF Regional Hospital &gt; Pediatrics</td><td>2</td><td>100%</td></tr><tr><td>South West &gt; California &gt; Sacramento &gt; Sacramento Regional Hospital &gt; Pediatrics</td><td>2</td><td>100%</td></tr><tr><td>South West &gt; California &gt; San Diego &gt; San Diego Regional Hospital &gt; Pediatrics</td><td>2</td><td>100%</td></tr></tbody></table>								Enterprise Unit Path	# with Version	% with Version	South West > California > Los Angeles > LA Regional Hospital > Pediatrics	3	100%	South West > California > San Francisco > SF Regional Hospital > Pediatrics	2	100%	South West > California > Sacramento > Sacramento Regional Hospital > Pediatrics	2	100%	South West > California > San Diego > San Diego Regional Hospital > Pediatrics	2	100%
Enterprise Unit Path	# with Version	% with Version																				
South West > California > Los Angeles > LA Regional Hospital > Pediatrics	3	100%																				
South West > California > San Francisco > SF Regional Hospital > Pediatrics	2	100%																				
South West > California > Sacramento > Sacramento Regional Hospital > Pediatrics	2	100%																				
South West > California > San Diego > San Diego Regional Hospital > Pediatrics	2	100%																				
<input type="checkbox"/>	Very Long Library Name	Evo IQ	Syringe	1.2.3.4.5.6.7.8.9.10	100000	John Jacob Smith	MM/DD/YYYY															

10 Items per page < 1 2 3 ... 8 9 10 > 400 Total items Next Step

© 2022 Baxter 3.0.1 | Privacy Policy | Terms & Conditions

## IQE Home Navy Blue


The first thing to notice is the design contains an expandable table row, while what has been built does not. We thought this would give an opportunity to explore additional ideas.

The updated design follows the white background for the headline, gray body section, white pagination and gray footer.

Updates:

- We restructured the configuration and search areas
- There is vertical spacing between table rows, this makes it easier to scan and looks “cleaner”

# Reports



Current IQE  
(cropped view)

# Reports

**Reports**

Compliance Limits **Dose Rate Change** Device Usage

Library: GatewayCQI Test Start Date: 01-Jan-2022 End Date: 05-Dec-2022

Quick Search Search By [ ] Apply Save/Manage Filters Clear Filters

Advanced Filters

Change Type Dose Rate Increases Sort By Dose Rate Increases

**All Selected Drugs**

78% infusions had at least one dose rate change 2,150 of 2,756 infusions

3,256 (100%) total dose rate changes over 2,150 infusions with dose rate changes.

979 (30%) Pullbacks 847 (26%) Overrides 1,432 (44%) Within Limits

**Dose Change Programming Mode - All Selected Drugs**

- Auto Program 1,658 71%
- Manual 690 29%
- DERS 1,869 80%
- Basic Mode 479 20%

<input type="checkbox"/>	Care Area	Programming Mode	Dose Changes
<input type="checkbox"/>	Med Surg	All	123
<input type="checkbox"/>	Pediatrics	All	157
<input type="checkbox"/>	Maternity	All	207
<input type="checkbox"/>	Cardiac Care	All	203

## IQE Baxter Blue

There were extensive design updates on these screens.

- The libraries and dates have edit icons to make things less cluttered
- All options are now “inline”
- All colors come from the branding guidelines
- The “All Selected Drugs” area is outlined with a white background to make it stand out more
- The Dose rate change flyout has a white background with a drop shadow
  - The checkboxes have rounded edges
  - Table no longer has zebra striping
  - Added sort functionality to table

# Reports

**Reports**

Compliance Limits **Dose Rate Change** Device Usage

Library: GatewayCQI Test Start Date: 01-Jan-2022 | End Date: 05-Dec-2022

Quick Search Search By  **Apply** [Save/Manage Filters](#) [Clear Filters](#)

Advanced Filters

Change Type Dose Rate Increases  Sort By Dose Rate Increases

**All Selected Drugs**

78% infusions had at least one dose rate change 2,150 of 2,756 infusions

3,256 (100%) total dose rate changes over 2,150 infusions with dose rate changes.

979 (30%) Pullbacks 847 (26%) Overrides 1,432 (44%) Within Limits

**Dose Change Programming Mode - All Selected Drugs**

- Auto Program**  
1,658 71%
- Manual**  
690 29%
- DERS**  
1,869 80%
- Basic Mode**  
479 20%

<input type="checkbox"/> Care Area	<input type="text"/> Programming Mode	<input type="text"/> Dose Changes
<input type="checkbox"/> Med Surg	All	<a href="#">123</a>
<input type="checkbox"/> Pediatrics	All	<a href="#">157</a>
<input type="checkbox"/> Maternity	All	<a href="#">207</a>
<input type="checkbox"/> Cardiac Care	All	<a href="#">203</a>


## IQE Home Navy Blue

There were extensive design updates on these screens.

- The libraries and dates have edit icons to make things less cluttered
- All options are now “inline”
- All colors come from the branding guidelines
- The “All Selected Drugs” area is outlined with a white background to make it stand out more
- The Dose rate change flyout has a white background with a drop shadow
  - The checkboxes have rounded edges
  - Table no longer has zebra striping
  - Added sort functionality to table

# Dark Mode Concepts

We have been noticing “dark mode” is becoming more prevalent. We decided it would make sense to attempt some dark mode options. These are based on a combination of the new Jellyfish interface, combined with the Fugu layouts using the new branding colors. We did several options using Baxter blue, Navy, gray and Cyan for buttons. We used Cyan for text links as it showed up the best. Here are a series of different variations...we went a little overboard on the options...





# Dark Mode Home

**Baxter** IQ Enterprise CONNECTIVITY SUITE JD

## IQ Enterprise Connectivity Suite

IQ Enterprise Connectivity Suite is a centralized platform for Baxter applications and devices that help inform, improve and deliver hospital therapies.

### Applications

ALL INFUSION PUMPS Search

- Configuration Management**  
Configure device distributions
- Device Management**  
Device registration & Management
- Reports**  
CQI Reports and Analytics
- Enterprise Hierarchy**  
Manage Enterprise Hierarchy

© 2022 Baxter 3.0.1 | [Privacy Policy](#) | [Terms & Conditions](#)

Baxter Blue

# Dark Mode Home

**Baxter**

**IQ Enterprise Connectivity Suite**  
IQ Enterprise Connectivity Suite is a centralized platform for Baxter applications and devices that help inform, improve and deliver hospital therapies.

**Applications** **ALL** INFUSION PUMPS Search

- Configuration Management**  
Configure device distributions
- Device Management**  
Device registration & Management
- Reports**  
CQI Reports and Analytics
- Enterprise Hierarchy**  
Manage Enterprise Hierarchy

© 2022 Baxter 3.0.1 | [Privacy Policy](#) | [Terms & Conditions](#)

Baxter Blue 75% opacity

# Dark Mode Home

The screenshot shows the dark mode interface of the IQ Enterprise Connectivity Suite. At the top left is the Baxter logo. The top right contains a help icon, the product name 'IQ Enterprise CONNECTIVITY SUITE', and a user profile icon 'JD'. Below the header is a title 'IQ Enterprise Connectivity Suite' and a descriptive sentence. A main section titled 'Applications' features a filter for 'INFUSION PUMPS' and a search bar. Four application cards are displayed: 'Configuration Management' (Configure device distributions), 'Device Management' (Device registration & Management), 'Reports' (CQI Reports and Analytics), and 'Enterprise Hierarchy' (Manage Enterprise Hierarchy). The footer includes copyright information and links to the privacy policy and terms & conditions.

**Baxter**

IQ Enterprise Connectivity Suite

IQ Enterprise Connectivity Suite is a centralized platform for Baxter applications and devices that help inform, improve and deliver hospital therapies.

**Applications** ALL INFUSION PUMPS Search

- Configuration Management**  
Configure device distributions
- Device Management**  
Device registration & Management
- Reports**  
CQI Reports and Analytics
- Enterprise Hierarchy**  
Manage Enterprise Hierarchy

© 2022 Baxter 3.0.1 | [Privacy Policy](#) | [Terms & Conditions](#)

Gray

# Dark Mode Home

**Baxter**

IQ Enterprise Connectivity Suite

IQ Enterprise Connectivity Suite is a centralized platform for Baxter applications and devices that help inform, improve and deliver hospital therapies.


**Applications** ALL INFUSION PUMPS Search

- Configuration Management**  
Configure device distributions
- Device Management**  
Device registration & Management
- Reports**  
CQI Reports and Analytics
- Enterprise Hierarchy**  
Manage Enterprise Hierarchy

© 2022 Baxter 3.0.1 | [Privacy Policy](#) | [Terms & Conditions](#)


Cyan

# Dark Mode Home: Navigation



White Links







# Dark Mode Home: Navigation



White Links

# Dark Mode Home: Navigation

## Variants



# Dark Mode Configuration Management

**Baxter** IQ Enterprise CONNECTIVITY SUITE JD

## Configuration Management

Configurations Distributions

Select Configuration Type: Drug Library Search By: Search By Search

Library Name	Device Model	Device Type	Library Version	Devices in Use	Imported By	Import Date
Very Long Library Name	Evo IQ	Syringe	1.2.3.4.5.6.7.8.9.10	100000	John Jacob Smith	MM/DD/YYYY
Very Long Library Name	Evo IQ	Syringe	1.2.3.4.5.6.7.8.9.10	100000	John Jacob Smith	MM/DD/YYYY

Enterprise Unit Path # with Version % with Version

- South West > California > Los Angeles > LA Regional Hospital > Pediatrics 3 100%
- South West > California > San Francisco > SF Regional Hospital > Pediatrics 2 100%
- South West > California > Sacramento > Sacramento Regional Hospital > Pediatrics 2 100%
- South West > California > San Diego > San Diego Regional Hospital > Pediatrics 2 100%

10 Items per page 1 2 3 ... 8 9 10 400 Total items Next Step

© 2022 Baxter 3.0.1 | Privacy Policy | Terms & Conditions

Baxter Blue 75% opacity buttons



# Dark Mode Configuration Management

**Baxter** IQ Enterprise CONNECTIVITY SUITE JD

## Configuration Management

Configurations Distributions

Select Configuration Type: Drug Library Search By: Search By Search

Library Name	Device Model	Device Type	Library Version	Devices in Use	Imported By	Import Date
Very Long Library Name	Evo IQ	Syringe	1.2.3.4.5.6.7.8.9.10	100000	John Jacob Smith	MM/DD/YYYY
Very Long Library Name	Evo IQ	Syringe	1.2.3.4.5.6.7.8.9.10	100000	John Jacob Smith	MM/DD/YYYY

Enterprise Unit Path

Enterprise Unit Path	# with Version	% with Version
South West > California > Los Angeles > LA Regional Hospital > Pediatrics	3	100%
South West > California > San Francisco > SF Regional Hospital > Pediatrics	2	100%
South West > California > Sacramento > Sacramento Regional Hospital > Pediatrics	2	100%
South West > California > San Diego > San Diego Regional Hospital > Pediatrics	2	100%

10 Items per page 1 2 3 ... 8 9 10 400 Total items Next Step

© 2022 Baxter 3.0.1 | Privacy Policy | Terms & Conditions

Baxter Blue buttons

# Dark Mode Configuration Management

The screenshot displays the Baxter Configuration Management interface in dark mode. At the top, the Baxter logo is on the left, and the user profile 'JD' and 'IQ Enterprise CONNECTIVITY SUITE' are on the right. The main heading is 'Configuration Management'. Below it, there are tabs for 'Configurations' and 'Distributions'. A search bar is present with a 'Search' button. The main content area shows a table of configurations with columns: Library Name, Device Model, Device Type, Library Version, Devices in Use, Imported By, and Import Date. A modal window is open, showing a table of 'Enterprise Unit Path' with columns for '# with Version' and '% with Version'. The footer contains copyright information, version '3.0.1', and links for 'Privacy Policy' and 'Terms & Conditions'.

Library Name	Device Model	Device Type	Library Version	Devices in Use	Imported By	Import Date
Very Long Library Name	Evo IQ	Syringe	1.2.3.4.5.6.7.8.9.10	100000	John Jacob Smith	MM/DD/YYYY
Very Long Library Name	Evo IQ	Syringe	1.2.3.4.5.6.7.8.9.10	100000	John Jacob Smith	MM/DD/YYYY

Enterprise Unit Path	# with Version	% with Version
South West > California > Los Angeles > LA Regional Hospital > Pediatrics	3	100%
South West > California > San Francisco > SF Regional Hospital > Pediatrics	2	100%
South West > California > Sacramento > Sacramento Regional Hospital > Pediatrics	2	100%
South West > California > San Diego > San Diego Regional Hospital > Pediatrics	2	100%

Gray buttons

# Dark Mode Configuration Management

**Baxter** IQ Enterprise CONNECTIVITY SUITE JD

## Configuration Management

Configurations Distributions

Select Configuration Type: Drug Library Search By: Search By Search

Library Name	Device Model	Device Type	Library Version	Devices in Use	Imported By	Import Date
Very Long Library Name	Evo IQ	Syringe	1.2.3.4.5.6.7.8.9.10	100000	John Jacob Smith	MM/DD/YYYY
Very Long Library Name	Evo IQ	Syringe	1.2.3.4.5.6.7.8.9.10	100000	John Jacob Smith	MM/DD/YYYY

Enterprise Unit Path # with Version % with Version

- South West > California > Los Angeles > LA Regional Hospital > **Pediatrics** 3 100%
- South West > California > San Francisco > SF Regional Hospital > **Pediatrics** 2 100%
- South West > California > Sacramento > Sacramento Regional Hospital > **Pediatrics** 2 100%
- South West > California > San Diego > San Diego Regional Hospital > **Pediatrics** 2 100%

10 Items per page 1 2 3 ... 8 9 10 400 Total items Next Step

© 2022 Baxter 3.0.1 | [Privacy Policy](#) | [Terms & Conditions](#)

Cyan buttons

# Reports

Compliance Limits **Dose Rate Change** Device Usage

Library: GatewayCQI Test Start Date: 01-Jan-2022 End Date: 05-Dec-2022

Quick Search Search By  **Apply** [Save/Manage Filters](#) [Clear Filters](#)

Advanced Filters


Change Type **Dose Rate Increases** Sort By **Dose Rate Increases**

### All Selected Drugs


78% infusions had at least one dose rate change 2,150 of 2,756 infusions

3,256 (100%) total dose rate changes over 2,150 infusions with dose rate changes.

979 (30%) Pullbacks 847 (26%) Overrides 1,432 (44%) Within Limits



### Dose Change Programming Mode - All Selected Drugs



<input type="checkbox"/>	Care Area	Programming Mode	Dose Changes
<input type="checkbox"/>	Med Surg	All	123
<input type="checkbox"/>	Pediatrics	All	157
<input type="checkbox"/>	Maternity	All	207
<input type="checkbox"/>	Cardiac Care	All	203
<input type="checkbox"/>	Care Area	All	157
<input type="checkbox"/>	Care Area	All	157
<input type="checkbox"/>	Care Area	All	157
<input type="checkbox"/>	Care Area	All	157
<input type="checkbox"/>	Care Area	All	157
<input type="checkbox"/>	Care Area	All	157

2,756 Total Infusions

**Export**

# Dark Mode Reports

Baxter Blue 75% opacity buttons

# Reports

Compliance Limits **Dose Rate Change** Device Usage

Library: GatewayCQI Test Start Date: 01-Jan-2022 End Date: 05-Dec-2022

Quick Search Search By  **Apply** [Save/Manage Filters](#) [Clear Filters](#)

Advanced Filters

Change Type **Dose Rate Increases** Sort By **Dose Rate Increases**

### All Selected Drugs


78% infusions had at least one dose rate change 2,150 of 2,756 infusions

3,256 (100%) total dose rate changes over 2,150 infusions with dose rate changes.

979 (30%) Pullbacks 847 (26%) Overrides 1,432 (44%) Within Limits



### Dose Change Programming Mode - All Selected Drugs



<input type="checkbox"/>	Care Area	Programming Mode	Dose Changes
<input type="checkbox"/>	Med Surg	All	123
<input type="checkbox"/>	Pediatrics	All	157
<input type="checkbox"/>	Maternity	All	207
<input type="checkbox"/>	Cardiac Care	All	203
<input type="checkbox"/>	Care Area	All	157
<input type="checkbox"/>	Care Area	All	157
<input type="checkbox"/>	Care Area	All	157
<input type="checkbox"/>	Care Area	All	157
<input type="checkbox"/>	Care Area	All	157
<input type="checkbox"/>	Care Area	All	157

2,756 Total Infusions

**Export**

# Dark Mode Reports

Baxter Blue buttons

No background on bar charts

# Reports

Compliance Limits **Dose Rate Change** Device Usage

Library: GatewayCQI Test Start Date: 01-Jan-2022 End Date: 05-Dec-2022

Quick Search Search By Apply Save/Manage Filters Clear Filters

Advanced Filters


Change Type Dose Rate Increases Sort By Dose Rate Increases

### All Selected Drugs


78% infusions had at least one dose rate change 2,150 of 2,756 infusions

3,256 (100%) total dose rate changes over 2,150 infusions with dose rate changes.

979 (30%) Pullbacks 847 (26%) Overrides 1,432 (44%) Within Limits



### Dose Change Programming Mode - All Selected Drugs



Care Area	Programming Mode	Dose Changes
<input type="checkbox"/> Med Surg	All	123
<input type="checkbox"/> Pediatrics	All	157
<input type="checkbox"/> Maternity	All	207
<input type="checkbox"/> Cardiac Care	All	203
<input type="checkbox"/> Care Area	All	157
<input type="checkbox"/> Care Area	All	157
<input type="checkbox"/> Care Area	All	157
<input type="checkbox"/> Care Area	All	157
<input type="checkbox"/> Care Area	All	157
<input type="checkbox"/> Care Area	All	157

2,756 Total Infusions

Export

# Dark Mode Reports

Gray buttons

# Reports

Compliance Limits Dose Rate Change Device Usage

Library: GatewayCQI Test Start Date: 01-Jan-2022 End Date: 05-Dec-2022

Quick Search Search By   [Save/Manage Filters](#) [Clear Filters](#)

Advanced Filters


Change Type Dose Rate Increases Sort By Dose Rate Increases

## All Selected Drugs


78% infusions had at least one dose rate change 2,150 of 2,756 infusions

3,256 (100%) total dose rate changes over 2,150 infusions with dose rate changes.

979 (30%) Pullbacks 847 (26%) Overrides 1,432 (44%) Within Limits



## Dose Change Programming Mode - All Selected Drugs



Care Area	Programming Mode	Dose Changes
<input type="checkbox"/> Med Surg	All	123
<input type="checkbox"/> Pediatrics	All	157
<input type="checkbox"/> Maternity	All	207
<input type="checkbox"/> Cardiac Care	All	203
<input type="checkbox"/> Care Area	All	157
<input type="checkbox"/> Care Area	All	157
<input type="checkbox"/> Care Area	All	157
<input type="checkbox"/> Care Area	All	157
<input type="checkbox"/> Care Area	All	157
<input type="checkbox"/> Care Area	All	157

2,756 Total Infusions

# Dark Mode Reports

Cyan buttons

No background, only outline on flyout

The word "Baxter" is written in a white, bold, italicized sans-serif font, centered on a background of overlapping blue and white geometric shapes.

***Baxter***

***Thank you!***